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- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih W.J, Bachorik P.S, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3351 – 364
- Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

HDL Cholesterol. Direct. Enzymatic colorimetric
Quantitative Determination of HDL Cholesterol
Only for *in vitro* use in clinical laboratory
Store at 2-8°C

Ref.: HDL-010/ HDL-110
HDL-010B / HDL-010S
1x60+1x20 mL/ 3x125+1x125 mL
1x750mL + 1x250 mL / 1x30+1x10 mL
CAL: 1x 1 mL

HDL CHOLESTEROL

PRINCIPLE OF THE METHOD

Directly determination of serum HDLc (high-density lipoprotein cholesterol) levels without the need for any pre-treatment or centrifugation of the sample^{3,5}.

The assay takes place in two steps.

– 1° Elimination of lipoprotein non-HDL

Cholesterol esters $\xrightarrow{\text{CHE}}$ Cholesterol + Fatty acids

Cholesterol + O₂ $\xrightarrow{\text{CHOD}}$ Cholestenone + H₂O₂

2 H₂O₂ $\xrightarrow{\text{Catalase}}$ 2 H₂O + O₂

– 2° Measurement of HDLc

Cholesterol esters $\xrightarrow{\text{CHE}}$ Cholesterol + Fatty acids

Cholesterol + O₂ $\xrightarrow{\text{CHOD}}$ Cholestenone + H₂O₂

2 H₂O₂ + HDAOS + 4-AA $\xrightarrow{\text{POD}}$ Quinone Pigment + 4H₂O

The intensity of the color formed is proportional to the HDLc concentration in the sample.

CLINICAL SIGNIFICANCE

HDL particles serve to transport lipoproteins in the blood-stream.

HDL is known as "good cholesterol" because high levels are thought to lower the risk of heart disease and coronary artery disease.

A low HDL cholesterol levels, is considered a greater heart disease risk^{1,2,4}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	N,N-bis(2-hydroxyethyl)-2-aminoethanesulphonic acid pH 6,6	100 mM
	N-(2-hydroxy-3-sulfofpropyl)-3,5-dimethoxyaniline (HDAOS)	0.7 mM
	Cholesterol Esterase	≥ 800 U/L
	Cholesterol oxidase	≥ 500 U/L
	Catalase	≥ 8300 KU/L
R 2	Ascorbic oxidase	≥ 3000 U/L
	N,N-bis(2-hydroxyethyl)-2-aminoethanesulphonic acid pH 7,0	1,1 mol/L
	4 – Aminoantipyrine (4-AA)	100 mM
	Peroxidase	≥ 30500 U/L
HDLc/ LDLc CAL		
Calibrator: Lyophilized human serum.		

PRECAUTIONS

HDLc/ LDLc CAL

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.

PREPARATION

– R 1 and R 2: Are ready to use.

– HDLc/ LDLc CAL: Dissolve the contents with 1 mL of distilled water. Cap vial and mix gently to dissolve contents.

STORAGE AND STABILITY (Note 1)

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not freeze the reagents.

– HDLc/ LDLc CAL: Once reconstitute 30 hours at 20-25°C, 2 weeks at 2-8°C or 3 months at -20°C.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

– Presence of particles and turbidity.

ADDITIONAL EQUIPMENT

– Spectrophotometer or colorimeter measuring at 600 nm.

– Matched cuvettes 1,0 cm light path.

– General laboratory equipment.

SAMPLES

Serum, heparinized plasma or EDTA plasma.

If any samples show precipitates, centrifuge before using.

Serum stable for 6 days at 2-8°C and 1 year when stored at -70°C.

PROCEDURE

- Assay conditions:
Wavelength: 600 -700 nm
Cuvette: 1 cm light path
Temperature: 37°C
- Adjust the instrument to zero with distilled water.
- Pipette into a cuvette:

	Blank	Calibrator	Sample
R 1 (μL)	300	300	300
Calibrator (μL)	--	3	--
Sample (μL)	--	--	3

- Mix and incubate for 5 min at 37°C and read the absorbance (A₁) of the samples and calibrator.

- Add:

	Blank	Calibrator	Sample
R 2 (μL)	100	100	100

- Mix and incubate for 5 min. at 37°C and read the absorbance (A₂) of the samples and calibrator, against the Blank.

CALCULATIONS

$\frac{(\Delta A) \text{ Sample} - (\Delta A) \text{ Blank}}{(\Delta A) \text{ Calibrator} - (\Delta A) \text{ Blank}} \times \text{Calibrator conc.} = \text{mg/dL of HDL-c in the sample}$

sample

Conversion factor: mg/dL x 0.0259= mmol/L.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: Control H Normal and Pathologic.

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES

	Men	Women
Low risk	> 50 mg/dL	> 60 mg/dL
Normal risk	35 – 50 mg/dL	45 – 60 mg/dL
High risk	< 35 mg/dL	< 45 mg/dL

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 5.0 mg/dL to linearity limit of 151 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (mg/dL)	28.0	76.1	27.5	75.3
SD	0.25	0.81	1.26	2.04
CV (%)	0.89	1.06	4.60	2.71

Sensitivity: 1mg/dL = 0,001399 (A).

Accuracy: Results obtained using BSM reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0.938.

Regression equation: y= 0.9825x – 1.41606 .

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

No interferences were observed to bilirubin, up to 30 mg/dL, hemoglobin up to 500 mg/dL, rheumatoid factors up to 1000 IU/mL or lipemia up to 1200 mg/dL.

Lipemic samples with a triglyceride concentration >1200 mg/dL should be diluted 1/10 with NaCl 9 g/L and multiply the result by 10.

NOTES

The reagent 2 presents yellowish coloration due to the peroxidase, but it does not affect its functionality.

BSM has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3351 – 364
- Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.