Bilirrubina Total. DMSO.

Ref.: BTC-021 /BTC-021B

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BIBLIOGRAFIA

- Kaplan A et al. Bilirubin. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1238-1241. 436 and 650.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2: 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966; Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Total Bilirubin. DMSO.

Ref.: BTC-021 /BTC-021B

Quantitative determination of total bilirubin total

Only for in vitro use in clinical laboratory

Store at 2-8°C

R1: 2x125 mL / 1x1000 mL

R2: 1x10 mL / 4x10 mL

CAL: 1x1 mL

TOTAL BILIRRUBIN**PRINCIPLE OF THE METHOD**

Bilirubin is converted to colored azobilirubin by diazotized sulfanilic acid and measured photometrically. Of the two fractions presents in serum, bilirubin-glucuronide and free bilirubin loosely bound to albumin, only the former reacts directly in aqueous solution (bilirubin direct), while free bilirubin requires solubilization with dimethylsulphoxide (DMSO) to react (bilirubin indirect). In the determination of indirect bilirubin the direct is also determined, the results correspond to total bilirubin.

The intensity of the color formed is proportional to the bilirubin concentration in the sample^{1,2,3}.

CLINICAL SIGNIFICANCE

Bilirubin is a breakdown product of hemoglobin, insoluble in water. It is transported from the spleen to the liver and excreted into bile.

Hyperbilirubinemia results from the increase of bilirubin concentrations in plasma.

Causes of hyperbilirubinemia:

Total bilirubin: Increase hemolysis, genetic errors, neonatal jaundice, ineffective erythropoiesis, and drugs.

Direct bilirubin: Hepatic cholestasis, genetic errors, hepatocellular damage^{1,6,7}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REACTIVOS

R 1	Sulfanilic acid	30 mmol/L
	Hydrochloric acid (HCl)	50 mmol/L
	Dimethylsulphoxide (DMSO)	7 mol/L
R 2	Sodium nitrite	29 mmol/L
Optional	BILIRUBIN CAL	

PRECAUTIONS

R1: H290-May be corrosive to metals. H314-Causes severe burns and eye damage. EUH208-Contains sulphanilic acid. May produce an allergic reaction.

Follow the precautionary statements given in MSDS and label of the product

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Color development in R 2.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 555 nm.

- Matched cuvettes 1,0 cm light path.

- General laboratory equipment.

SAMPLES

Serum or plasma, free of hemolysis¹. Protect samples from direct light.

Stability: Bilirubin is stable at 2-8°C for 4 days and 2 months at -20°C.

PROCEDURE

- Assay conditions:
Wavelength: 555 nm (530-580)
Cuvette: 1 cm light path
Temperature: 15-25°C
- Adjust the instrument to zero with distilled water.
- Pipette into a cuvette:
-

	Blank	B. Total
R 1 (mL)	1,5	1,5
R 2 (μL)	--	50
Sample ^(Note 1) / Calibrator (μL)	100	100

- Mix and incubate for exactly **5 minutes** at room temperature.
- Read the absorbance (A).

CALCULATIONS**a. With Calibrator:**

$$\frac{(A) \text{ Sample} - (A) \text{ Sample Blank}}{(A) \text{ Calibrator} - (A) \text{ Calibrator Blank}} \times \text{Conc. Calibrator} = \text{mg/dL}$$

b. With Factor:

$$(A) \text{ Sample} - (A) \text{ Sample Blank} \times \text{Factor}^* = \text{mg/dL}$$

$$\text{*Factor: } \frac{\text{Concentration of Calibrator}}{(A) \text{ Calibrator} - (A) \text{ Calibrator Blank}}$$

Conversion factor: mg/dL x 17.1 = μmol/L.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures. If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹

Up to 1,10 mg/dL $\equiv$ Up to 18,81 μmol/L

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 0,00526 mg/L to *linearity limit* of 18 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (mg/dL)	1,53	5,06	1,53	5,02
SD	0,03	0,05	0,03	0,11
CV (%)	1,73	1,01	1,92	2,18

Sensitivity: 1 mg/dL = 0.05074 A.

Accuracy: Results obtained using BSM reagents (y) did not show systematic differences when compared with other commercial reagents (x). The results obtained using 50 samples were the following:
Correlation coefficient (r)²: 0.991.

Regression equation: y = 0,82743x + 0,0382

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

Hemolysis causes decreased bilirubin values^{1,2,3}.

A list of drugs and other interfering substances with bilirubin has been reported by Young et. al^{4,5}.

NOTES

- For bilirubin determination in newborns, pipette 50 μL of sample. Multiply the result by 2.
- BSM has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

- Kaplan A et al. Bilirubin. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1238-1241. 436 and 650.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966; Acta 13; 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

Total Bilirubin. DMSO.

Quantitative determination of total bilirubin total
Only for in vitro use in clinical laboratory
Store at 2-8°C

Ref.: BTC-021 /BTC-021B

R1: 2x125 mL / 1x1000 mL
R2: 1x10 mL / 4x10 mL
CAL: 1x1 mL

TOTAL BILIRRUBIN